

# SINGLE END SCALERS

Single end scalers are an excellent choice for pedodontic practices to remove deposits from supragingival surfaces and in orthodontic practices to remove cement. Single end scalers should be used when a double end instrument is not suitable or required. Nordent will produce any of our hygiene tip patterns in a single end configuration upon request.



## Sickle N5

A moderate size, medium-reach sickle with a blade width of 0.9mm that tapers to a point.

Handle Selection:	<input checked="" type="radio"/> CESC N5 (shown)	<input type="radio"/> RESC N5	<input type="radio"/> SC N5
-------------------	--	-------------------------------	-----------------------------



## Sickle 109

A moderate size, long-reach sickle with a blade width of 0.9 mm that tapers to a point.

Handle Selection:
<input checked="" type="radio"/> CESC109
<input type="radio"/> RESC109
<input type="radio"/> SC109



## Towner U15

A large size, long-reach sickle with a blade width of 1 mm that tapers to a point.

Handle Selection:
<input checked="" type="radio"/> CESC U15
<input type="radio"/> RESC U15
<input type="radio"/> SC U15



## Ivory CI-1

A moderate size, medium-reach sickle with a blade width of 2 mm that tapers to a point.

Handle Selection:
<input checked="" type="radio"/> CESC CI1
<input type="radio"/> RESC CI1
<input type="radio"/> SC CI1



## Younger-Good 15

A small size, short-reach sickle on a straight shank with a blade width of 1 mm that tapers to a point.

Handle Selection:
<input checked="" type="radio"/> CESC YG15
<input type="radio"/> RESC YG15
<input type="radio"/> SC YG15



## Univ. of So. Cal. 128

A moderate size, medium-reach sickle on a straight shank with a blade width of 1.2 mm that tapers to a point.

Handle Selection:
<input checked="" type="radio"/> CESC 128
<input type="radio"/> RESC 128
<input type="radio"/> SC 128



## Whiteside 2

A small size, short-reach sickle on a straight shank with a 1.5 mm blade that tapers to a point.

Handle Selection:
<input checked="" type="radio"/> CESC W2
<input type="radio"/> RESC W2
<input type="radio"/> SC W2



## Nebraska 128

Thinner version of the USC128 with a blade width of 1mm.

Handle Selection:
<input checked="" type="radio"/> CESC NE128
<input type="radio"/> RESC NE128
<input type="radio"/> SC NE128



## Jacquette J1S

A small size, straight blade Jacquette with a 0.9 mm blade width that tapers to a point. This tip is also known as the Jacquette #33.

Handle Selection:
<input checked="" type="radio"/> CESC J1S
<input type="radio"/> RESC J1S
<input type="radio"/> SC J1S



## "B" Scaler

A moderate size, straight blade with a 0.5 mm width that has two sharp edges on the face and the back of the blade. This enables the blade to be used in a push or pull stroke.

Handle Selection:
<input checked="" type="radio"/> CESC B
<input type="radio"/> RESC B
<input type="radio"/> SC B