

GOLDMAN-FOX

The Goldman-Fox series of instruments includes universal curettes, Jacquettes, sickles and hoes.

Scalers & Curettes



GF #1 Anterior Jacquette The blades are 5 mm long and taper to a point from a 1 mm width.

Handle Selection: CESC GF1 (shown) RESC GF1 ESC GF1 RSC GF1 SC GF1



GF #2 Anterior Curette Designed like a Gracey medium contra-angle with an open blade curve and a rigid shank. The blade is 1.1 mm wide with a round toe.

Handle Selection: CESC GF2 (shown) RESC GF2 ESC GF2 RSC GF2 SC GF2



GF #3 Universal Curette The terminal shank is 1.1 mm in diameter and 10 mm long. The blade width is 1 mm.

Handle Selection: CESC GF3 (shown) RESC GF3 ESC GF3 RSC GF3 SC GF3



GF #4 Universal Curette The terminal shank is 1.1 mm in diameter and 12 mm long. The blade width is 1 mm.

Handle Selection: CESC GF4 (shown) RESC GF4 ESC GF4 RSC GF4 SC GF4



GF #5 Hoe Mirror image hoes that are 1.5 mm wide. For facial and lingual surfaces of posterior teeth.

Handle Selection: CESC GF5 (shown) RESC GF5 ESC GF5 RSC GF5 SC GF5



GF #6 Hoe The blades are 1.5 mm wide and are well-suited for anterior areas.

Handle Selection: CESC GF6 (shown) RESC GF6 ESC GF6 RSC GF6 SC GF6



GF #21 Anterior Scaler Combines a small Jacquette that has a 3.5 mm length and a small sickle that has a 1.2 mm blade that tapers to a point.

Handle Selection: CESC GF21 (shown) RESC GF21 ESC GF21 RSC GF21 SC GF21