

# ORTHODONTIC INSTRUMENTS



**Band Pusher – Black's #6-7** Two identical oval-shaped, serrated blades set at opposing angles. The blades measure 3.5 mm wide and 1.5 mm thick.

Handle Selection: ● CNB6-7



**Band Pusher – Scaler #114**

Handle Selection: ◻ EOTN114



**Scaler – Ligature Director #115**

Handle Selection: ◻ EOTN115



**Band Pusher – Scaler #119**

Handle Selection: ◻ EOTN119



**Ligature Director #120 Single End**

Handle Selection: ● OTN120

**Mathieu #207C**

Medium jaws with carbide inserts and fine serrations. The unique ratchet mechanism allows the instrument to be locked and opened by simply squeezing the handles. Excellent for left-handed operation (5"/130 mm).

