



MICRO SURGERY

Tunnel Surgery Instruments

The Nordent Tunnel Surgery Instruments are designed by Prof. Dr. Ates Parlar, Professor of Periodontology, Ankara University and also in private practice in Ankara, Turkey.



The instruments feature very thin tapered blades, oval in shape, that are sharp around their entire periphery. Blade detail is shown at 300% of actual size.



Papilla Retractor #1 Parlar Papilla retractor universal, blue color rings.

Handle Selection:  CEEPAP1 (shown)



Papilla Retractor #2 Parlar Papilla retractor buccal offset 45°, white color rings.

Handle Selection:  CEEPAP2 (shown)



Papilla Retractor #3 Parlar Papilla retractor lingual offset 45°, red color rings.

Handle Selection:  CEEPAP3 (shown)



Papilla Retractor #4 Parlar Papilla retractor sharp periosteal, lingual/buccal blades, green color rings.

Handle Selection:  CEEPAP4 (shown)

Castroviejo Scissors #6536

Castroviejo scissors 18 cm. Curved blades, round European style handle.

S6536 (shown at 75% actual size)



Castroviejo Needle Holder #5333

Castroviejo needle holder 18 cm. Carbide beaks, round European style handle.

NH5333 (shown at 75% actual size)

