

NEEDLE HOLDERS

Castroviejo

Castroviejo needle holders feature beaks with very fine serrations on narrow profile jaws. The locking mechanism is activated by squeezing the spring-action handle. All edges of the hinges are rounded off and tapered to eliminate any snagging of the suture material.

Shown at 50% actual size unless noted.

Stainless Steel Jaws

Carbide Jaws



actual size



Castroviejo #209

Very fine serrations on narrow profile jaws. The locking mechanism is activated by squeezing the spring-action handle. This is a great needle holder for periodontal applications where a very fine suture is used (5" / 130 mm).

NH209 (actual size)



actual size

Castroviejo #208

Very fine serrations on narrow profile jaws with carbide inserts. The locking mechanism is activated by squeezing the spring-action handle. This is a great needle holder for periodontal applications where a very fine suture is used (5" / 130 mm).

NH208



actual size

Castroviejo NH4062

Very fine serrations on narrow profile jaws with carbide inserts. Curved beaks (5.5" / 140 mm).

NH4062



actual size

Castroviejo NH4388

Very fine serrations on narrow profile jaws with carbide inserts. Straight beaks, round Euro style handle (5.5" / 140 mm).

NH4388



actual size

Castroviejo NH4389

Very fine serrations on narrow profile jaws with carbide inserts. Curved beaks, round Euro style handle (5.5" / 140 mm).

NH4389

NEEDLE HOLDERS

Nordent stainless steel and carbide needle holders feature fine machined serrations on the beaks. The edges of the hinges are tapered and rounded off to prevent any snagging of the suture material. The carbide needle holders can be identified by their gold handles and are recommended for use for those practices with a high volume of surgical procedures.

Shown at 50% actual size unless noted.

Carbide Jaws




Mayo-Hegar #213
Large jaws with carbide inserts and fine serrations (7" / 180 mm).

NH213 (shown at 75% actual size)



Mayo-Hegar #201
Large jaws with carbide inserts and fine serrations (6 1/2" / 160 mm).

NH201



Crile-Wood #203
Medium jaws with carbide inserts and fine serrations (6" / 150 mm).

NH203



Derf #205
Medium jaws with carbide inserts and fine serrations (5" / 130 mm).

NH205

Stainless Steel Jaws




Mayo-Hegar #202
Large jaws with fine serrations (6 1/2" / 160 mm).

NH202



Crile-Wood #204
Medium jaws with fine serrations (6" / 150 mm).

NH204



Derf #206
Medium jaws with fine serrations (5" / 130 mm).

NH206