

# RETRACTORS

Used for retracting cheek and tongue during surgical procedures. All stainless steel and can be sterilized by any method.



**University of Minnesota**  
(5 3/4" / 145 mm)

*RCM (shown at 50% actual size)*



**Weider #28**  
28 mm wide "Sweetheart" retractor.  
(5 1/4" / 135 mm)

*RCW28 (shown at 50% actual size)*



**Weider #36**  
36 mm wide "Sweetheart" retractor.  
(5 3/4" / 145 mm)


*RCW36 (shown at 50% actual size)*

# SUCTION TIPS

Nordent Suction Tips are all stainless steel and are used to evacuate fluids from the surgical site. Frazier Suction Tips have a relief hole on the thumb rest that allows the practitioner to regulate the intensity of the suction at the tip. Each comes with its own stylette that is used to clear debris from the inside of the tube.

**Frazier #10**  
2.0 mm internal diameter

*ATF10 (shown at 75% actual size)*



**Frazier #12**  
3.0 mm internal diameter

*ATF12 (shown at 75% actual size)*

