

RONGEURS

Used to contour and remove alveolar bone during extraction procedures, Nordent Rongeurs are hardened and tempered to maintain their sharpness and have a satin finish to eliminate glare.



Friedman #1

Narrow profile blades cut around the entire periphery and are set at a 30° angle (5-1/2"/140 mm).

R1



Mini-Friedman #1A

Narrow profile blades cut around the entire periphery and are set at a 30° angle (5" / 125mm).

R1A



Blumenthal #2

Standard profile blades cut around the entire periphery and are set at a 30° angle (6-1/4" / 160 mm).

R2



Cleveland #4

Standard profile blades are set at a 30° angle (5-1/2" / 140 mm).

R4

Did you know...

Nordent will sharpen and repair all brands of Rongeurs. Call us today for more information
1-800-966-7336