

ROOT TIP PICKS

Root picks are used to retrieve loose root fragments after extraction. Nordent uses a special heat-treating process that assures a sharp, long-lasting point and shank strength. The handle is all stainless steel and can be sterilized by any method.

Double-ended root tip picks and root teasers simplify tray setups for more efficient procedures. Nordent offers these instruments in large diameter Duralite® ColorRings™ handle or in the standard handle.



Heidbrink #1
Straight.

Heidbrink #2
Left.

Heidbrink #3
Right.

Heidbrink #2-3
Double-ended root pick with pointed tips.

Davis #11
Root tip teaser with a rounded tip.

Davis #11P
Root tip teaser with a pointed tip and serrated blade.

<p>ER1 (shown) ER1SR (serrated)</p>	<p>ER2 (shown) ER2SR (serrated)</p>	<p>ER3 (shown) ER3SR (serrated)</p>	<p>☉ CEER2-3 (shown) ● ER2-3</p>	<p>☉ CEER11 (shown) ● ER11</p>	<p>☉ CEER11P (shown) ● ER11P</p>
---	---	---	--------------------------------------	------------------------------------	--------------------------------------