

PERIODONTAL

Periodontal Knives

Nordent periodontal knives are forged from high-carbon stainless steel and hardened to the highest degree. Each blade is honed by expert technicians for a long-lasting, ultra-sharp edge.



Goldman Fox #11 Spear-shaped interdental knife with offset blades

Handle Selection: ● KGF11



Kirkland #15-16 Pointed kidney-shaped blades for gingivectomy bevel incisions.

Handle Selection: ● KK15-16



Orban #1-2 Double-ended mirror image blades on a contra-angled shank that are pointed for access to interproximal tissue.

Handle Selection: ● KOR1-2



Orban #1 Single-end (left).

Handle Selection: ● KOR1



Orban #2 Single-end (right).

Handle Selection: ● KOR2

Scalpel Blade Handles

Nordent scalpel blade handles are precision machined so that scalpel blades are held firmly in place. All handles are stainless steel so they can be sterilized by any method.



Bard Parker type with metric ruler (5"/125 mm).

HB3



Perfectly balanced round handle allows easy rotation in difficult-to-reach areas (6"/155 mm).

HB5