

SCISSORS

Shown at 50% actual size.

Iris



Iris #301
Straight blades
(4 1/2" / 115 mm).

S301



Iris #302
Curved blades
(4 1/2" / 115 mm).

S302

Castroviejo



Castroviejo #311
Curved blades with spring action
handle. Excellent for tissue and
suture (4" / 100 mm).

S311

LaGrange



LaGrange #314
A unique curvature allows easy
access to the tightest areas. One
blade is serrated to hold tissue
when cutting (4 1/2" / 115 mm).

S314

Suture



Suture #304
Angular blades with suture hook
(4 1/2" / 115 mm).

S304



Spencer Suture #306
Straight blades with suture hook
(4 3/4" / 120 mm).

S306



Suture #322
Long handle, straight blades with
suture hook (6 1/4" / 160 mm).

S322

Sullivan



Sullivan #316
Extra-long curved blades. One
blade is serrated to hold tissue
when cutting (5 1/2" / 140 mm).

S316

SCISSORS

Shown at 50% actual size.

Goldman-Fox



Goldman Fox #313
Curved blades with one blade serrated to hold tissue when cutting (5" / 130 mm).

S313



Goldman Fox #321
A unique handle curvature allows easy access to the tightest areas. One blade is serrated to hold tissue when cutting (5" / 130 mm).

S321

Kelly



Kelly #317
Long handle with straight blades (6 1/4" / 160 mm).

S317



Kelly #318
Long handle with curved blades (6 1/4" / 160 mm).

S318



Kelly #319
Long handle with angular blades (6 1/4" / 160 mm).

S319

Dean



Dean Gum #320
Long curved handle with straight blades that are angled upward. One blade is serrated to hold tissue when cutting (6 3/4" / 170 mm).

S320

Super-Cut

All Nordent super-cut scissors feature one blade that is honed thin with a knife edge. The opposing blade is a standard scissors blade with the addition of fine serrations. One finger ring is gold plated for easy identification.



Super-Cut - Iris Curved
Iris scissors, straight blades (4 1/2" / 115 mm).

S302SC



Super-Cut - LaGrange
Sharply curved blades and S curve handle (4 1/2" / 115 mm).

S314SC



Super-Cut - Goldman Fox
Curved tapered blades (5" / 130 mm).

S313SC