

# NORDENT® ONPOINT SCALERS AND CURETTES

## Featuring Nordent UltraSteel™ Technology

**BONUS:** Each instrument comes with a packet of Nordent ColorRings™ to customize your instrument to your own personalization

### Affordable

Receive the longest-lasting scalers and curettes in the dental market for an unbeatable price.


### Sharper

Features Nordent's exclusive UltraSteel™ technology. The hardest, sharpest, longest-lasting scalers and curettes available.

### Customizable

Nordent OnPoint handles come with two interchangeable color rings, one on each end of the instrument.

**Keep your Hygiene Practice  
OnPoint All Day Everyday**



# OnPoint Scalars and Curettes

## Nordent UltraSteel™ Technology

UltraSteel technology takes advantage of advancements in metallurgy and cryogenic hardening to produce a micro-structure not previously possible. All other scalars are made from stainless steels that have large, isolated carbide molecules. This is why hardness and wear resistance can be inconsistent.

UltraSteel technology features extremely fine carbide particles that are evenly distributed throughout the micro-structure to provide you with a blade that is consistently harder, sharper and longer lasting than ANY other brand.



UltraSteel Technology  
Scalars & Curettes

Under **500x magnification** these SEM photos clearly show the **superior micro-structure** and **uniform carbide particle distribution** of our **Nordent UltraSteel technology**



Leading sharpen-less scalar



All other major brands (440A)

## AVAILABLE NORDENT ONPOINT SCALARS AND CURETTES



### UNIVERSAL CURETTES - PURPLE

**Scalette N135**  
Posterior - Universal OPN135

**Barnhart #5-6**  
Universal OPBH5-6

**McCall #13S-14S**  
Universal OPMC13S-14S

**Columbia #13-14**  
Universal OPC013-14



### GRACEY CURETTES - TEAL

**Gracey 11-12**  
Posterior OPGR11-12

**Gracey 13-14**  
Posterior OPGR13-14



### POSTERIOR SCALARS - LIGHT PINK

**Sickle #204S**  
Posterior OP204S

**Sickle #204SD**  
Posterior OP204SD

**Offset Universal #2**  
Anterior OPN2



### UNIVERSAL SCALARS - YELLOW

**NMJ Scaler**  
(Universal R138) -  
Anterior OPNMJ



### ANTERIOR SCALARS - FUSCHIA

**Sickle / Spoon N5-N1**  
Anterior OPN5-N1

**Offset Sickle N67**  
Anterior OPN67

**Anterior Lingual  
Stain Remover**  
Anterior OPALSR



A passion for premium craftsmanship and industry-leading innovations has transformed Nordent over the last 50 years into the market leader you know today. Behind the sharpest name in dentistry is a team of trusted advisors dedicated to your success and ready to provide you with unbeatable support.



**PREMIUM  
CRAFTSMANSHIP**



**TRUSTED  
ADVISOR**



**UNBEATABLE  
SUPPORT**

© 2023 Nordent and all other marks are registered trademarks of Nordent Manufacturing Inc. All rights reserved. Product and specifications are subject to change without notice.



T: 800.966.7336 US & Canada | T: 847.437.4780 | F: 847.437.4786  
customerservice@nordent.com | www.nordent.com

