

CUSTOM CASSETTES

Maximize efficiency
while protecting your
instruments

Increase your productivity,
saving up to six minutes
per procedure

Heavy gauge
stainless steel



Snap-In-Silicone™ Rails

Keep your kits organized and secure
with Snap-in-Silicone rails

- Protect sterilized instruments
- Reduce risk of sharps injuries
- Incur less breakage

FEATURES




Custom Cassettes


- Heavy gauge stainless steel doesn't flex or bend
- Customizable to your instrument kit
- Unique hole design for optimal flow
- Electro-polished finish for stain resistance
- Save processing time with double-sided cassettes

Easy Crescent Glide Latch

- Fast and easy to operate, open and go
- Secure and durable



Easy to Open
and Close



Innovative Snap-In-Silicone™ Rails

- Uniquely designed to easily secure instruments
- Optimized instrument spacing for fast, one-handed retrieval and return
- Increase instrument exposure during sterilization with minimal contact
- Available in seven colors for easy organization

**5 YEAR UNCONDITIONAL
GUARANTEE**

Flexible Designs, Premium Craftsmanship

Custom Cassettes with Snap-In-Silicone™ rails can be tailored to your instrument kit. Listed below are the most popular configurations.

Small Cassette 3" W X 8" L X 1.25" H

Available configurations:
0-8 instruments

8 INSTRUMENTS

- CD1-8N1 Blue
- CD1-8N2 Teal
- CD1-8N3 Fuchsia
- CD1-8N4 Yellow
- CD1-8N5 Light Green
- CD1-8N6 White
- CD1-8N7 Purple



CD1-8N1 Shown

Medium Cassette 5.5" W X 8" L X 1.25" H

Available configurations:
0-16 instruments

8 INSTRUMENTS

- CD2-8D1 Blue
- CD2-8D1M Blue (Mat)
- CD2-8D2 Teal
- CD2-8D3 Fuchsia
- CD2-8D4 Yellow
- CD2-8D5 Light Green
- CD2-8D6 White
- CD2-8D7 Purple

16 INSTRUMENTS

- CD2-16N1 Blue
- CD2-16N2 Teal
- CD2-16N3 Fuchsia
- CD2-16N4 Yellow
- CD2-16N5 Light Green
- CD2-16N6 White
- CD2-16N7 Purple



CD2-16N1 Shown

Large Cassette 8" W X 11" L X 1.25" H

Available configurations:
0-32 instruments

16 INSTRUMENTS

- CD3-16D1 Blue
- CD3-16D1M Blue (Mat)
- CD3-16D2 Teal
- CD3-16D3 Fuchsia
- CD3-16D4 Yellow
- CD3-16D5 Light Green
- CD3-16D6 White
- CD3-16D7 Purple

24 INSTRUMENTS

- CD3-24D1 Blue
- CD3-24D1M Blue (Mat)
- CD3-24D2 Teal
- CD3-24D3 Fuchsia
- CD3-24D4 Yellow
- CD3-24D5 Light Green
- CD3-24D6 White
- CD3-24D7 Purple



CD3-16D1 (Mat) Shown



Single and additional double-sided configurations are also available

Learn more at www.nordent.com or by contacting us at **800-966-7336**

ACCESSORIES


Snap-In-Silicone™ Rails

Available in 8 (4 top, 4 bottom) and 16 (8 top and 8 bottom). Colors available include: blue, teal, fuchsia, yellow, white, light green, and purple.




Silicone Finger Mats

Available in small, medium and large mats for larger or bulky items. Mats can be cut to your exact requirements. Only available in blue.




Air Water Syringe Clips

Secure your air and water syringe tips, scissors and needle holder to the inside of the cassette.



Large and Small Parts Boxes

Keep small items safe and secure during sterilization.



Contact one of your **Nordent trusted advisors** today to learn how to incorporate cassettes based on your specific needs.



A passion for premium craftsmanship and industry-leading innovations has transformed Nordent over the last 50 years into the market leader you know today. Behind the sharpest name in dentistry is a team of trusted advisors dedicated to your success and ready to provide you with unbeatable support.



**PREMIUM
CRAFTSMANSHIP**



**TRUSTED
ADVISOR**



**UNBEATABLE
SUPPORT**

© 2022 Nordent and all other marks are registered trademarks of Nordent Manufacturing Inc. All rights reserved. Nordent Custom Cassettes US Patent No. 11,324,577



T: 800.966.7336 US & Canada | T: 847.437.4780 | F: 847.437.4786
customerservice@nordent.com | www.nordent.com

