

# UltraSteel™ Technology




## The hardest, longest-lasting edge possible<sup>1</sup>

### IT'S WHAT'S ON THE INSIDE THAT COUNTS!

Nordent's exclusive UltraSteel technology takes advantage of advancements in metallurgy and cryogenic hardening to produce a micro-structure not previously possible. All other scalers are made from stainless steels that have large, isolated carbide molecules. This is why hardness and wear resistance can be inconsistent. UltraSteel technology features extremely fine carbide particles that are evenly distributed throughout the micro-structure to provide you with a blade that is consistently harder, sharper and longer lasting than **ANY** other brand.


Under 500x magnification these SEM photos clearly show the superior micro-structure and uniform carbide particle distribution of our Nordent UltraSteel technology



XDURA® & OnPoint™  
Scalers and Currettes



Leading sharpen-less scaler



All other brands (440A)

### Featured in XDURA Instruments



### Featured in OnPoint Instruments



To learn more about XDURA and OnPoint Scalers and Currettes, visit [Nordent.com](http://Nordent.com)

© 2024 Nordent and all other marks are registered trademarks of Nordent Manufacturing Inc. All rights reserved. Product and specifications are subject to change without notice.

1. Data on file by request



T: 800.966.7336 US & Canada | T: 847.437.4780 | F: 847.437.4786

