LUXPOINT ULTRASONIC INSERTS



Featuring Nordent UltraSteel™ Technology

Elevate your scaling with Nordent® LuxPoint™ ultrasonic inserts that feature our exclusive UltraSteel™ technology, providing enhanced performance and industry-leading tip durability.



Proven Durability

UltraSteel technology is known for its consistent hardness ensuring the long lasting durability of LuxPoint tips. Thanks to UltraSteel superior microstructure and uniform carbide particle distribution, LuxPoint tips last over 50% longer than other leading brands.

Reliable Performance

Our universal tips fit any ultrasonic scaler unit, providing seamless performance.

Comfort Meets Precision

Our unique hexagon grip pattern offers enhanced ergonomics to minimize strain and maximize productivity. The silicon grip enhances your sense of touch, giving you comfort and a secure hold.

Enhanced Patient Experience

Our direct tip targeted water flow limits unnecessary spray, offering better visibility and a more comfortable experience for patients.

*Data on File by Request

FEATURES AT A GLANCE











90-day Warranty

LUXPOINT™ - Ultrasonic Inserts

MODELS AVAILABLE



#10 Universal Slim 30kHz

Color: Nordent Light Purple

SKU: UIUS3010



#10 Universal Slim 25kHz

Color: Nordent Light Purple

SKU: UIUS2510



#1000 Triple Bend Slim 30kHz

Color: Nordent Dark Purple

SKU: UIUS301000



#1000 Triple Bend Slim 25kHz

Color: Nordent Dark Purple

SKU: UIUS251000



#10 Universal X-Slim 30kHz

Color: Nordent Blue

SKU: UIUS30XT

Universal Slim

An all-around great choice for

Universal Extra Thin

tip design improves access.

Triple Bend Slim

triple bend design an excellent choice access. Also great for orthodontic





A passion for premium craftsmanship and industry-leading innovations has transformed Nordent over the last 55 years into the market leader you know today. Behind the sharpest name in dentistry is a team of trusted advisors dedicated to your success and ready to provide you with unbeatable support.





© 2024 Nordent and all other marks are registered trademarks of Nordent Manufacturing Inc. All rights reserved. Product and specifications are subject to change without notice.





