


# LUXPOINT™ ULTRASONIC INSERTS

**Nordent®**  
The sharpest name in dentistry.

Featuring Nordent UltraSteel™ Technology

Elevate your scaling with Nordent® LuxPoint™ ultrasonic inserts that feature our exclusive UltraSteel™ technology, providing enhanced performance and industry-leading tip durability.



Lasts over **50% longer** than other industry leading ultrasonic inserts\*

## Proven Durability

UltraSteel technology is known for its consistent hardness ensuring the long lasting durability of LuxPoint tips. Thanks to UltraSteel superior microstructure and uniform carbide particle distribution, LuxPoint tips last over 50% longer than other leading brands.

## Reliable Performance

Our universal tips fit any ultrasonic scaler unit, providing seamless performance.

## Comfort Meets Precision

Our unique hexagon grip pattern offers enhanced ergonomics to minimize strain and maximize productivity. The silicon grip enhances your sense of touch, giving you comfort and a secure hold.

## Enhanced Patient Experience

Our direct tip targeted water flow limits unnecessary spray, offering better visibility and a more comfortable experience for patients.

\*Data on File by Request

## FEATURES AT A GLANCE



Direct Tip  
Water Flow



Easy-to-identify  
Grip Colors



Precision Control  
and Comfort




Compatible  
Performance



90-day  
Warranty

# LUXPOINT™ - Ultrasonic Inserts

## MODELS AVAILABLE



**#10 Universal Slim 30kHz**

**Color:** Nordent Light Purple


**SKU:** 146-3611



**#10 Universal Slim 25kHz**

**Color:** Nordent Light Purple


**SKU:** 146-3612



**#1000 Triple Bend Slim 30kHz**

**Color:** Nordent Dark Purple


**SKU:** 146-3613



**#1000 Triple Bend Slim 25kHz**

**Color:** Nordent Dark Purple

**SKU:** 146-3614



**#10 Universal X-Slim 30kHz**

**Color:** Nordent Blue

**SKU:** 146-3615

### Universal Slim

An all-around great choice for a variety of scenarios, in a slim design for maximum subgingival and interproximal access.

### Universal Extra Thin

Designed for subgingival and interproximal surfaces with light to moderate deposits. The thin tip design improves access. The recommended power level is low-medium.

### Triple Bend Slim

Recommended for light to heavy deposits. Beveled edges make the triple bend design an excellent choice for stubborn calculus. The triple bend shank allows increased access around line angles and interproximal access. Also great for orthodontic cement removal. Recommended for low-high power or medium high for tenacious calculus.



A passion for premium craftsmanship and industry-leading innovations has transformed Nordent over the last 55 years into the market leader you know today. Behind the sharpest name in dentistry is a team of trusted advisors dedicated to your success and ready to provide you with unbeatable support.



**PREMIUM CRAFTSMANSHIP**



**TRUSTED ADVISOR**



**UNBEATABLE SUPPORT**

\*Data on File by Request

© 2024 Nordent and all other marks are registered trademarks of Nordent Manufacturing Inc. All rights reserved. Product and specifications are subject to change without notice.

T: 800.966.7336 US & Canada | T: 847.437.4780 | F: 847.437.4786  
customerservice@nordent.com | www.nordent.com

