

HARTZELL™ CONE SOCKET

SCALERS AND CURETTES



Signature stainless steel or sharpen-free titanium nitride (TiN) coated tips

Lightweight, ergonomically designed 9.5 mm diameter handle

Secure, comfortable handle grip

The original Hartzell™ cone socket instrument that you know and trust is back!

Save Costs


Replace one tip at a time when needed with Hartzell cone socket instruments, saving on the lifetime cost of your scalers and curettes.

Environmentally Friendly

Disposing of the tips instead of the full instrument allows your dental office to reduce waste compared to using conventional instruments.

Mix & Match

Customize your Hartzell handle with your preferred instrument tips: Our tried-and-true signature stainless steel tips or long-lasting, sharpen-free titanium nitride (TiN) infused tips. Both are available to suit your unique preferences.




The original Hartzell designs are now masterfully produced by the craftspeople at Nordent®

Available Hartzell™ Cone Socket Scalers and Curettes



TIP DESCRIPTION	440A	TiN
Hartzell UC/Ratcliff 1, C	R1S	R1ST
Hartzell UC/Ratcliff 2, C	R2S	R2ST
Hartzell Ratcliff 3, C	R3S	R3ST
Hartzell Ratcliff 4, C	R4S	R4ST
Hartzell Anterior Offset Sickle N6, S	S6S	S6ST
Hartzell Anterior Offset Sickle N7, S	S7S	S7ST
Hartzell Universal Offset Sickle NMJ Left, S	MJ1S	MJ1ST
Hartzell Universal Offset Sickle NMJ Right, S	MJ2S	MJ2ST
Hartzell Gracey 5 Rigid, C	GL5S	GL5ST
Hartzell Gracey 6 Rigid, C	GL6S	GL6ST
Hartzell Gracey 7 Rigid, C	GL7S	GL7ST
Hartzell Gracey 8 Rigid, C	GL8S	GL8ST
Hartzell Gracey 11 Rigid, C	GL11S	GL11ST
Hartzell Gracey 12 Rigid, C	GL12S	GL12ST
Hartzell Gracey 13 Rigid, C	GL13S	GL13ST
Hartzell Gracey 14 Rigid, C	GL14S	GL14ST
Hartzell Columbia 4 Right, C	4RS	4RST
Hartzell Columbia 4 Left, C	4LS	4LST
Hartzell Posterior Sickle 204S Right, S	204RS	204RST
Hartzell Posterior Sickle 204S Left, S	204LS	204LST
Hartzell Cone Socket Handle	CSHA	
Hartzell Wrench	WRHZ	



“I’m thrilled to partner with the Nordent Team to reintroduce the iconic Hartzell designs to the dental industry.”

–Andrew McIver-Hartzell



****All Hartzell Tips Attach to CSHA Handle**

C= Curette

S= Scaler



A passion for premium craftsmanship and industry-leading innovations has transformed Nordent over the last 55 years into the market leader you know today. Behind the sharpest name in dentistry is a team of trusted advisors dedicated to your success and ready to provide you with unbeatable support.



PREMIUM CRAFTSMANSHIP



TRUSTED ADVISOR



UNBEATABLE SUPPORT

© 2024 Nordent and all other marks are registered trademarks of Nordent Manufacturing Inc. All rights reserved. Product and specifications are subject to change without notice.



T: 800.966.7336 US & Canada | T: 847.437.4780 | F: 847.437.4786
customerservice@nordent.com | www.nordent.com