

# TITANIUM IMPLANT PROBES



## GENTLE & LONG-LASTING TIPS

Lightweight DuraLite® Round handle is made with stainless steel for better tactile sensitivity and comfort.



### Color-Coded Tips

Color-coding is applied to Probes under extreme temperatures so that the color penetrates the surface of the tip, preventing the markings from easily scratching off.



Our titanium tips are an ideal alternative to plastic or stainless steel, offering an implant-friendly option.

#### 1. **REPCT33**

#N33 (Color-Coded) Probe, SE,  
DuraLite Round Handle (3-6-9-12 mm)

#### 2. **REPCT15**

UNC #N15 (Color-Coded) Probe, SE,  
DuraLite Round Handle (1-15 mm)