

ONPOINT[®]

SCALERS AND CURETTES

PREMIUM DESIGN, PRICED FOR YOUR PRACTICE



Nordent **ColorRings[®]** provide you with a fun way to personalize your hygiene instruments.

Discover **enhanced comfort** with the **10 mm updated handle design**.

Nordent **UltraSteel[™]** technology for a long-lasting edge.

Premium Quality

Nordent tips made in the USA by master craftsmen.


Sharper

UltraSteel[™] technology results in tips with long-lasting sharpness and strength.

Customizable

Nordent OnPoint handles come with two interchangeable ColorRings[®], one on each end of the instrument.

OnPoint[®]



Explore our wide range of patterns, including our extensive selection of Gracey patterns!

ONPOINT® SCALERS AND CURETTES

Nordent UltraSteel™ Technology


UltraSteel technology takes advantage of advancements in metallurgy and cryogenic hardening to produce a micro-structure not previously possible. All other scalers are made from stainless steels that have large, isolated carbide molecules. This is why hardness and wear resistance can be inconsistent.

UltraSteel technology features extremely fine carbide particles that are evenly distributed throughout the micro-structure, resulting in tips with long-lasting sharpness and strength.


Under **500x magnification**, these SEM photos clearly show the **superior micro-structure** and **uniform carbide particle distribution** of our Nordent UltraSteel technology.



UltraSteel Technology
Scalers & Curettes



Leading sharpen-less scaler



All other major brands (440A)

Extensive Selection of Gracey Instruments

Whether your preference is Nordent's own TRUGracey® or traditional Gracey patterns, our selection meets all your patient needs.




OPPGR1-2



OPPGR7-8



OPPGR11-12



OPPGR13-14

View our full section of scaler and curette patterns at nordent.com



A passion for premium craftsmanship and industry-leading innovations has transformed Nordent over the last 55 years into the market leader you know today. Behind the sharpest name in dentistry is a team of trusted advisors dedicated to your success and ready to provide you with unbeatable support.



**PREMIUM
CRAFTSMANSHIP**



**TRUSTED
ADVISOR**



**UNBEATABLE
SUPPORT**

© 2025 Nordent and all other marks are either registered or owned trademarks of Nordent Manufacturing, Inc. All rights reserved. Product and specifications are subject to change without notice. DTD #OP7125