

GRACEY & TRUGRACEY®

CURETTES

Choose what works best for your style and your patients.

Nordent Gracey Patterns

Nordent Gracey curettes feature a more common blade design and are widely known as “Area Specific Curettes.” The straighter blades allow for a mix of various strokes geared towards everyday hygiene.

Nordent TRUGracey® Patterns

Nordent’s exclusive TRUGracey curettes are designed like Dr. Gracey intended them to be. Our blade design conforms to the true intent of Dr. Gracey, with blades that wrap around the tooth anatomy and easily slide into the gingival pocket with a vertical stroke.



COMPARE FOR YOURSELF

| Features | Nordent Gracey Curettes | Nordent TRUGracey® Curettes |
|--------------|---|---|
| Stroke Type | Clinicians using a mix of vertical, oblique, and horizontal strokes | Clinicians using primarily vertical strokes |
| Blade Shape | Standard blade with a gentle radius | True radius blade, more curved |
| Pocket Entry | Traditional entry with control across stroke types | Due to less overhang of the curette toe, enter using vertical technique |

Nordent Gracey




CESCGR1-2
(1 shown)

Nordent TRUGracey®



OPPGA1-2
(1 shown)



NORDENT GRACEY DESIGNS

Standard Gracey




Features a standard length terminal shank and blade.

Rigid Gracey




Features a rigid shank and thicker blade that are well-suited to periodontal and surgical application.

Long Reach Gracey



Features an extended terminal shank that is 3 mm longer than our Standard Gracey curettes.

Mini-Long Gracey



Features a shorter, narrower blade and an extended terminal shank that is 3 mm longer than Standard Gracey curettes.

NORDENT TRUGRACEY® DESIGNS

Standard TRUGracey



Features a blade face curve with a true radius that conforms to tooth anatomy better.

Long Reach TRUGracey



Features an extended terminal shank that is 3 mm longer than Standard TRUGracey.

| | XDURA® DuraLite® ColorRings® | XDURA® DuraLite® Round | OnPoint® | Relyant® | Hartzell™ | DuraLite® Round | DuraLite® ColorRings® | DuraLite® HEX | Medium Round |
|-------------------|------------------------------|------------------------|----------|----------|-----------|-----------------|-----------------------|---------------|--------------|
| Gracey | | | | | | | | | |
| Standard | ● | ● | ● | ● | | ● | ● | ● | ● |
| Rigid | | | | ● | ● | ● | ● | | ● |
| Long Reach | ● | ● | ● | ● | | ● | ● | ● | ● |
| Mini-Long | ● | ● | ● | ● | | ● | ● | ● | ● |
| TRUGracey® | | | | | | | | | |
| Standard | ● | ● | ● | ● | | ● | ● | ● | ● |
| Long Reach | ● | ● | ● | ● | | ● | ● | ● | ● |

Scan to View
Full Product Page



Complimentary sharpening
available for all Nordent XDURA®
and Relyant® instruments