


LUXPOINT™

ULTRASONIC INSERTS



Featuring Nordent UltraSteel™ Technology

Elevate your scaling with Nordent® LuxPoint™ ultrasonic inserts that feature our exclusive UltraSteel™ technology, providing enhanced performance and industry-leading tip durability.



Proven Durability

UltraSteel technology known for its consistent hardness ensures the extended durability of LuxPoint tips. Thanks to UltraSteel superior microstructure and uniform carbide particle distribution, LuxPoint tips result in long-lasting.

Reliable Performance

Our universal tips fit any ultrasonic scaler unit, providing seamless performance.

Comfort Meets Precision

Our unique hexagon grip pattern, offers enhanced ergonomics to minimize strain and maximize productivity. The silicon grip enhances your sense of touch, giving you comfort and a secure hold.

Enhanced Patient Experience

Our direct tip targeted water flow limits unnecessary spray, offering better visibility and a more comfortable experience for patients.

*Data on File by Request

FEATURES AT A GLANCE



Direct Tip
Water Flow



Easy-to-identify
Grip Colors



Precision Control
and Comfort



Compatible
Performance



90-day
Warranty

LUXPOINT™ - Ultrasonic Inserts

MODELS AVAILABLE



#10 Universal Slim 30kHz

Color: Nordent Light Purple


SKU: UIUS3010



#10 Universal Slim 25kHz

Color: Nordent Light Purple


SKU: UIUS2510



#1000 Triple Bend Slim 30kHz

Color: Nordent Dark Purple


SKU: UIUS301000



#1000 Triple Bend Slim 25kHz

Color: Nordent Dark Purple

SKU: UIUS251000



#10 Universal X-Slim 30kHz

Color: Nordent Blue

SKU: UIUS30XT

Universal Slim

An all-around great choice for a variety of scenarios, in a slim design for maximum subgingival and interproximal access.

Universal Extra Thin

Designed for subgingival and interproximal surfaces with light to moderate deposits. The thin tip design improves access. The recommended power level is low-medium.

Triple Bend Slim

Recommended for light to heavy deposits. Beveled edges make the triple bend design an excellent choice for stubborn calculus. The triple bend shank allows increased access around line angles and interproximal access. Also great for orthodontic cement removal. Recommended for low-high power or medium high for tenacious calculus.



A passion for premium craftsmanship and industry-leading innovations has transformed Nordent over the last 55 years into the market leader you know today. Behind the sharpest name in dentistry is a team of trusted advisors dedicated to your success and ready to provide you with unbeatable support.



PREMIUM CRAFTSMANSHIP



TRUSTED ADVISOR



UNBEATABLE SUPPORT

*Data on File by Request

© 2026 Nordent and all other marks are registered trademarks of Nordent Manufacturing Inc. All rights reserved. Product and specifications are subject to change without notice.

T: 800.966.7336 US & Canada | T: 847.437.4780 | F: 847.437.4786
customerservice@nordent.com | www.nordent.com

